



## Analysis of Beer

H. Daidoji,  
HORIBA Ltd, Japan

**Keywords:** food, beverage, method development  
**Elements:** Al, Ca, Cd, Cu, Fe, Mg, Mn, Na, Zn

### 1 Introduction

The analysis of beer is important to the brewery industry, both for production purposes as well as to the final consumer. The analysis can pose several problems due to both the CO<sub>2</sub> gas dissolved in the beer and the dissolved and undissolved protein. This report demonstrates several approaches and compares the results obtained.

### 2 Principle

#### 2.1 Technique used

The elemental analysis of these samples was undertaken by Inductively Coupled Plasma Optical Emission Spectrometry (ICP-OES). The sample is nebulized then transferred to an argon plasma. It is decomposed, atomized and ionized whereby the atoms and ions are excited. We measure the intensity of the light emitted when the atoms or ions return to lower levels of energy. Each element emits light at characteristic wavelengths and these lines can be used for quantitative analysis after a calibration.

#### 2.2 Wavelength choice

The choice of the wavelength in a given matrix can be made using the profile function, or by using "Win-Image", which is rapid semi-quantitative analysis mode using multiple wavelengths. The principle is the same in either case: record the scans of analytes at low concentration, and of the matrix. By superimposing the spectra, we see possible interferences.

### 3 Sample Preparation

The analysis of beer is affected by CO<sub>2</sub> and organic materials such as protein. The CO<sub>2</sub> must

be degassed, otherwise CO<sub>2</sub> breaks out of solution during the pumping stage and nebulization. There are several methods to achieve this including bubbling N<sub>2</sub> into the solution, ultrasonic vibration and heating.

We tried three methods for the analysis.

- A) Typical calibration curve
- B) Internal standard method
- C) Standard addition method

#### A) Typical calibration curve

To each 20 g or 40 g sample of beer, about 5 ml of concentrated HNO<sub>3</sub> was added to dissolve organic materials. It was heated at 80 °C for 3 minutes. The action of HNO<sub>3</sub> and heat is to degas the CO<sub>2</sub>. Distilled water was added to a total volume of 50 mL.

The standards for the calibration, for all elements, were made from Spex Certiprep Mix standard solutions.

#### B) Internal standard method

A vessel containing 20 g of beer was kept in an ultrasonic vibration bath (80W) for 1 hour to degas the CO<sub>2</sub> in the beer. Due to the alcohol content the specific gravity can change, so it is important to weigh the sample, rather than pipette a volume. Then, 5 ppm Yttrium was added as an internal standard element in this solution and distilled water was added to make a total weight of 50 g. The standards for calibration curves, for each element, was made from Spex Certiprep Mix standard solutions and 5 ppm Yttrium was added at the same concentrations as the samples. An alternative to adding the Y to the samples directly is to introduce 5ppm using a



"Y" piece. This ensures all solutions have identical concentration of internal standard.

### C) Standard addition method

Five vessels containing 20 g of beer were kept in an ultrasonic vibration bath (80 W) for 1 hour. Due to the alcohol content the specific gravity can change, so it is important to weigh the sample, rather than pipette a volume. After degassing of CO<sub>2</sub>, mixed standard solutions were added in these solutions to prepare 0, 0.02, 0.05, 0.1 and 10 ppm spikes of each element. Background corrections were done for all elements.

## 4 Instrument specification

The work was done on a JY ULTIMA and is also applicable to a JY ULTIMA 2 ICP spectrometer. The specifications are listed below.

**Table 1: Specification of spectrometer**

Parameters	Specifications
Mounting	Czerny-Turner
Focal length	1 m
Thermoregulation	Yes
Variable resolution	Yes
Nitrogen purge	Yes
Grating number of grooves	2.400 gr/mm
Orders	2
1st order resolution	0.005 nm
2nd order resolution	0.010 nm

**Table 2: Specification of RF Generator**

Parameters	Specifications
Type of generator	Solid state
Observation	Radial
Frequency	40.68 MHz
Control of gas flowrate	By computer
Control of pump flow	By computer
Cooling	Air

The operating conditions of the spectrometer are listed in Table 3.

**Table 3: Operating conditions**

Parameter	Condition
Generator power	1000 W
Plasma gas flowrate	12 L/min
Auxiliary gas flowrate	0.4 L/min
Sheath gas flowrate	0.25 L/min
Nebulizer	1 L/min
Sample uptake	1 mL/min
Type of nebulizer	Tangential (Burgener)
Type of spray chamber	Cyclonic
Argon humidifier	No
Injector tube diameter	3.0 mm

## 5 Wavelength selection and analytical conditions

**Table 4: Analytical conditions**

Mode of calculation	Entrance Slits (μm)	Exit Slits	Integration time (sec)	Replicates
Maximum Mode	20 μm	15 μm	3 s	3

## 6 Results

Profiles for each element in beer dissolved with HNO<sub>3</sub> (first method) are shown in Figure 1.

Calibration curves for Mn for the three methods are given Figure 2.

The results are shown in Table 4 for the three types of calibration.



Table 4: Results for three types of calibration

Element	Wavelength nm	A (ppm)		B (ppm)	C (ppm)
		a	b		
Al	396.1	0.080	0.082	0.078	0.072
Ca	422.6	28.8	29.4	29.6	27.2
Cd	228.8	< 0.0001	< 0.0001	< 0.0001	< 0.0001
Cu	324.7	0.042	0.040	0.042	0.050
Fe	259.9	0.056	0.051	0.059	0.066
Mg	285.2	48.1	52.2	53.5	53.7
Mn	257.6	0.075	0.074	0.076	0.071
Na	588.9	29.9	31.0	30.7	30.5
Zn	213.9	0.0052	0.0035	0.0050	0.0058

## 7 Conclusion

To analyze beer, it was shown in this Application Note that the three types of methods used to eliminate the CO<sub>2</sub> and protein gave equivalent results.









Figure 1: Profiles of elements in beer.  
Dissolution with HNO<sub>3</sub>






### Typical calibration curve

s: 0.014988

BEC: 938 ng/l

Correlation: 0.9998




### Internal standard method

s: 0

BEC: 49.7 ng/l

Correlation: 1



### Standard addition method

s:0

BEC: 15.4 µg/l

Correlation: 1




Figure 2: The calibration curves Mn 257.610 nm



**In the USA:**  
 Jobin Yvon Inc.  
 3880 Park Avenue  
 Edison, NJ 08820  
 Tel: 1-732-494-8660  
 Fax: 1-732-494-8796  
 E-mail:  
 emission@jyhoriba.com

**In France:**  
 Jobin Yvon S.A.S.  
 16-18, rue du Canal  
 91165 Longjumeau Cedex  
 Tel: (33) 1/64 54 13 00  
 Fax: (33) 1/69 09 90 88

**In Japan:**  
 Horiba Ltd.  
 2 Miyano Higashi, Kisshoin  
 Minami-ku, Kyoto 601-8510  
 TEL: (81) 75 313 8121  
 FAX: (81) 75 321 5725  
 www.horiba.com

**Germany:** (49) 89/46 23 17-0  
**Italy:** (39) 2/57 60 56 90  
**U.K.:** (44) 20/82 04 81 42

**China:** (86) 10/6880 0240  
**Spain:** (34) 91/724 16 57  
**Other Countries:** Contact JY S.A.S.

1-877-JYHORIBA

www.jyhoriba.com



HORIBA GROUP