

Excellence in Thin Film Metrology - 180 Years of Optical Expertise

Spectroscopic Ellipsometry



Explore the future


HORIBA GROUP

EMISSION • FLUORESCENCE • FORENSICS • GRATINGS & DEM • OPTICAL SPECTROSCOPY • RAMAN • THIN FILM




Film Thickness

High sensitivity to ultra-thin films



Change in signal for 10 Å thick monolayer at the Brewster angle




Optically pumped semiconductor laser structure

- Accurate thin film measurement from a few angstroms to several microns
- For single layer or complex multilayer stacks


— Ultra-thin monolayer
 - - - Fused silica substrate

High resolution monochromator configuration



Ultra-thick dielectric layer (32µm)

Complex multilayer stack analysis




Spectroscopic
 Ellipsometry :
 Extensive
 Capabilities



Optical Properties


- Refractive index (n) and extinction coefficient (k) from the far-UV to near-IR for complex materials, graded and anisotropic layers

(n,k) Evolution for a graded layer



ZnO thin film

Anisotropic layer characterization





Organic film


Advanced
 Metrology to

- Composition / crystallinity
- Microstructure
- Film uniformity by area and depth


Inhomogeneity over depth: material proportion determination using EMA




Ge concentration determination



Excellent correlation between SIMS and ellipsometry



Automated sample mapping for area uniformity



Strengthen your Competitive Edge



A Concentration of Advantages

- Non-destructive technique
- Highly accurate & reproducible
- Can be used in any transparent and semi-transparent medium
- No reference material necessary
- Very sensitive, especially to ultra-thin films (< 10 nm)
- Simultaneous multiple parameter determination
- Measures data at wavelength of interest




Chemistry & Biology

- Polymers
- Langmuir Blodgett films
- Liquid
- Lipids, protein adhesion



OLED application





NPB electroluminescent organic film on ITO on glass substrate



Optical Coatings

- Wave guides
- Laser mirrors, AR
- Photovoltaic applications




Anti-reflective coating application




Flat Panel Displays

- Oxides, nitrides (ITO, MgO, TaO_x, Al₂O₃, SiN_x...)
- Amorphous silicon, polysilicon
- Colour filters
- Resist






Flat Panel Display application



Dielectric Materials

- Nitrides, oxides, fluorides, carbides




Widest
Ellipsometry
Applications
Database

Three
Applications
Laboratories


Semiconductor Technology



- Oxides, nitrides, oxinitrides
- Thin NO capacitors, ONO, OPO, ONOPO
- SOI, SIMOX
- Metals (Ti, TiN, TaN, MoSi, WSi_x, Al, Cu, Co...)
- Photoresists, ARC
- High k, low k



ONO application




Alloys




- Si_(1-x)Ge_(x), Hg_(1-x)Cd_(x)Te
- Al_(x)Ga_(1-x)As, InGaAsP
- Superlattices and MQW

Determination of alloy gradient composition



III-V application




Real Time Process Control



- Thickness monitoring
- Growth and etch rates
- Endpoint detection
- Alloy composition
- Crystallinity
- Surface damage
- Contamination

Process control multilayer

Filter application



And Any Application That Needs Thin Film Characterization...



- Data storage
CDR, CDR-W, DVD-R, DVD-RW
Magneto-optic materials (Faraday or Kerr effect)
- Non linear optical devices
- Telecommunications...

for
Providing
Expert
Solutions



UVISEL - Bench top configuration

Flexibility Combined with High Performance

The UVISEL Spectroscopic Phase Modulated Ellipsometer is a unique instrument that delivers the highest accuracy and precision for demanding research and industrial QC applications.

It covers a wide spectral range from 190 to 1700 nm, with a complete application database included in the software.



UVISEL
Fully integrated configuration

• Ex-Situ Configuration

The UVISEL instrument is highly featured and integrates high stability electronic systems with advanced software capabilities.

It is possible to automate the instrument to enhance its performance and range of application by adding many available options.


When integrated into a cabinet the UVISEL addresses the needs of industrial research and process development. Its advantages are a small footprint, an easy installation and clean-room compatibility.

• In-Situ Configuration

By mounting the ellipsometer onto a process chamber the speed and stability of the in-situ UVISEL allows real-time monitoring and control of thin film deposition or etch processes with monolayer resolution.

The intrinsically high acquisition rate of the UVISEL allows sampling at millisecond time resolution, making it the ideal solution for in-situ monitoring and real-time control.

The in-situ UVISEL is routinely used in plasma deposition/etching, thermal oxidation, surface cleaning, implantation, corrosion, MBE, CVD, PVD, electrochemistry...



UVISEL coupled to CVD chamber



From Research to Industry

The Ultimate Tool for Powerful Research

SPECIFICATION TABLE			UVISEL				
● standard ● option			Uvisel	Uvisel FUV	Uvisel NIR	Uvisel ER	Uvisel MWL
Spectral Range							
Visible - 240 - 830 nm			x				x
FUV - 190 - 830 nm				x			x
NIR - 260 - 1700 nm					x		
ER - 190 - 1700 nm						x	
Configuration							
Ex-situ	Bench top	●	Typical table layout : 150 wide x 100 deep x 85 high in cm				
	Cabinet	●	Overall dimensions : 100 wide x 80 deep x 150 high in cm				
In-situ							●
Mechanical & Optical Parts							
Light Source	75 W Xe-lamp	●			●		●
	150 W Xe-lamp			●		●	●
Sample Stage	manual	●	Manual height (20 mm), tilt, theta adjustment				
	XY automatic	●	8" supplied as standard with others sizes upon request				
Goniometer	manual	●	Manually adjustable angle from 55° to 75° by step of 5°				
	automatic	●	Automatically adjustable angle from 40° to 90° by step of 0.01°				
Microspot	manual	●	3 positions : 0.05 - 0.1 - 1 mm				
	automatic	●	4 positions : 0.06 - 0.12 - 0.25 - 1.2 mm (up to 0.03 mm upon request)				

Thin Film Production Control

Fast, Accurate and Stable to Ensure High Yields in Quality and Quantity

The UT-300 and FF-1000 instruments have been developed to provide specific process control solutions for the semiconductor and flat panel display industries. These accurate, automated thin film metrology tools deliver both unique performance and proven reliability for on-line quality control of production processes.

Equipped with achromatic microspot optics, wafer handling system, autofocus and pattern recognition software the UT-300 - Fully Automatic Ultra Thin Film Analyzer - accurately characterizes demanding thin film structures with a throughput in excess of 100 wafers/hour. The instrument has a deep-UV option (190 nm) and is compatible with 6", 8" and 12" wafers.



UT-300

The FF-1000 features all of the advantages of the UVISEL and UT-300 with a fully automated large area sample stage able to accept samples up to 1000 mm x 1000 mm.

Two detection modes are available:

- scanning mode for highly accurate characterization of material physical properties (n, k, d, composition...)
- multiwavelength mode for fast and precise production control


Powerful software based on Windows™ allows automated measurement and analysis of single or complex multiple layer stacks.



FF-1000

On-Line Control Oriented Systems

Phase Modulated Ellipsometry



Unequaled Capabilities for Accurate Ultra-Thin Film Characterization

Ellipsometry is based on the measurement of the light polarization change upon reflection from a sample surface or interface.

The experimental data are usually expressed as two parameters Ψ and Δ , which are related to the Fresnel reflection coefficients by :

$$\rho = \frac{r_p}{r_s} = \tan \Psi e^{i\Delta}$$

These two coefficients contain information related to material optical properties and physical dimensions. Spectroscopic ellipsometry measures this complex ratio ρ as a function of wavelength.

Three Key Factors for Success

• The Most Accurate Measurement of Δ Parameter

Spectroscopic Phase Modulated Ellipsometers (SPME) use photoelastic devices to perform the polarization modulation without any mechanical movement, resulting in :

- Excellent signal/noise ratio from FUV to NIR
- No insensitive regions

This technology is consequently the best suited for accurate ultra-thin film measurement on transparent substrates.

Uniqueness of Spectroscopic Phase Modulated Ellipsometry:


Accurate measurement of Δ around 0° and 180°



Measurement on ultra-clean silica at the Brewster angle


• High stability

State-of-the-art control of PEM allows very high system stability.



• The Highest Data Acquisition Speed

The 50 kHz modulation frequency of the Photoelastic Modulator (PEM) allows ultra fast acquisition at up to 1 ms per point for applications such as dynamic studies and liquid - surface measurements. The very high speed of the UVISEL in combination with digital signal averaging provides significant advantages over conventional ellipsometers.




Commitment
to
Excellence

Performance
Specification

Advanced Spectroscopic Ellipsometry Software

Powerful DeltaPsi Windows™ based software makes full use of all the benefits provided by state-of-the-art JY ellipsometry hardware. The largest variety of advanced modelling functions gives research engineers the full performance of ellipsometric analysis. A simple user interface allows the operator to perform routine tasks very easily.



Variable angle measurements

Key Features

- Acquisition and analysis of ellipsometric, kinetic, transmission and reflection data
- Advanced mathematical fitting algorithms
- Bibliographic reference database - extendable
- Data and graphs easily transferred to Windows™ applications
- Import/export package functions for high flexibility in file manipulation

Enhanced Modelling Functions


- Graded layer
- Roughness or interface
- Composition / crystallinity
- Anisotropic layer
- Thickness uniformity
- Depolarization factor
- Complete library of material properties based on dispersion relations
- Automatic backside correction for thick transparent substrates
- Periodic structure
- BLMC (proprietary mathematical algorithm) for ultra-thin film applications
- Multiguess, multistart, multimodel, correlation, ...



ITO graded layer model taking into account backside reflections

Simple Automatic Operation to Meet Production Needs

- Recipe procedure: data acquisition, analysis and mapping routines
 - Fitting procedure
 - Pattern recognition function
 - 2D and 3D display views



Pattern recognition function

Defined and Optimized



Company History

Jobin Yvon SAS was created in 1819 and over the years distinguished scientists such as Fresnel, Babinet, Laurent, Fabry, Perot and Fizeau have worked in collaboration with our organization. Jobin Yvon provides state-of-the-art spectroscopic instrumentation and components and is today the leading supplier of optical spectrometer assemblies and subassemblies to original equipment manufacturers (OEMs). Our instruments respond to the needs of research and engineers wherever there is a need for accurate and reliable metrology tools.

The Thin Film Division was created by combining the cumulative optics experience of Jobin Yvon with the expertise and almost twenty years experience of SOFIE Instruments in thin film processing and plasma diagnostics.

By a process of continuous collaboration with major semiconductor manufacturers, research establishments and plasma processing equipment suppliers the Thin Film Division has developed a wide range of instruments dedicated to advanced thin film metrology such as Ellipsometry, Optical Emission Spectroscopy (OES) and Interferometry.

Our instruments are manufactured under a strict quality assurance program with the aim to satisfy customer needs by the highest level of quality and performance.

Pursuing a tradition of innovation Jobin Yvon acquired the license for a Phase Modulated Ellipsometer developed in collaboration with CNRS laboratories (license Anvar-Ecole Polytechnique).

In 1992 this collaboration won an Award for research-to-industry technology transfer, and Jobin Yvon successfully developed a broad variety of spectroscopic ellipsometry instruments from this base.

Jobin Yvon Spectroscopic Ellipsometers provide the highest level of performance currently available with continuous improvement in terms of hardware performance and analysis software.

This commitment to product excellence and continued support is part of JY's culture.



**Building
Relationships
with
Creative
Solutions**

**Reliability
and Continued
Support**

**We are committed to
give you the right answer
the first time.**

Call us and let us work together!

Service and Application Support

A Worldwide Experienced Service Team for Your Technical Support

Jobin Yvon offers a complete range of service and preventative maintenance plans to fit your needs and for your complete satisfaction.

Purposes of this contract are :

- System installation
- Technical training for a smooth and carefree start up
- Regular preventative maintenance actions to keep your instrument in its best working condition

Our staff of highly trained service and application engineers stand ready to provide assistance when and where you need it.

- Application training programs - advanced or basic -
- Training session seminars on site
- Sample analysis support for new applications
- Software configuration according to your applications
- Supplies, accessories and upgrades for the future growth of your instrument

Sample Measurement Service

Jobin Yvon offers sample measurement services in its application lab.

The laboratory is equipped with ellipsometers from FUV to NIR including all options.

An expert application team will characterize your samples and provide you with a detailed analysis report.



Our extensive experience in ellipsometry applications has allowed us to compile the widest analysis database.

Feel free to contact the sales division for a quotation, in relation with our application engineers for the most appropriate technical evaluation.

Global Network

To give you easy access to our technologies and services, a global network through subsidiaries and distributors operates in more than 80 different countries as well as three application laboratories, based in France, USA and Japan, for your sample analysis.

- Jobin Yvon Offices



**are Part
of JY's
Culture**

A Full Spectrum of Solutions

Emission

Spark, ICP and GDS spectrometers,
Particle size analyzer,
Gas analysis.

Fluorescence

Routine, modular spectrofluorometers,
Spectrofluorometers for lifetime and steady state
measurements.

Forensics

Forensic light sources.
Fingerprint instrumentation for detection and
comparison.

Gratings - OEM - VUV

Holographic and blazed holographic gratings,
OEM monochromator and spectrograph
module, Contract manufacturing.

VUV Instrumentation

Synchrotron beamlines,
Spectrographs Monochromators.

Thin Film

Advanced Film Thickness Metrology
Spectroscopic Ellipsometry.

Advanced End Point Detection

OES. Laser and white light interferometry.
Reflectometry.

Optical Spectroscopy

Wide range of monochromators (0.1 to 3 m focal length).
Light sources, spectrographs, CCD arrays accessories and
software.

Raman

Laboratory, research and process instruments from
UV to NIR.

In the USA:

Jobin Yvon Inc.
3880 Park Avenue
Edison, NJ 08820-3012
Toll-Free: **1-866-jobinyvon**
Tel: +1-732-494-8660
Fax: +1-732-549-5125
E-mail: info@jobinyvon.com
www.jobinyvon.com

In France:

Jobin Yvon S.A.S.
Z.A. de la Vigne aux Loups - 5, Av. Arago - 91380 Chilly Mazarin
Tel : +33 (0) 1 64 54 13 00 - Fax: +33 (0) 1 69 74 88 61
Email: tfd-sales@jobinyvon.fr
www.jobinyvon.fr

Germany: +49 (0) 89 4623 17-0

Italy: +39 0 2 57603050

U.K.: +44 (0) 20 8204 8142

In Japan:

Horiba Ltd., JY Optical Sales Dept.
Higashi-kanda Daiji Building
1-7-8 Higashi-kanda, chiyoda-ku
Tokyo 101-0031
Tel: +81 (0) 3 3861 8231
www.jyhoriba.jp

JY China: +86 (0)10 6849 2216

JY Korea: +82 2 576 8650



This document is not contractually binding under any circumstances - Printed in France - ADECOM +33 (0)1 64 53 11 22 - 02/2004

Explore the future

HORIBA GROUP

EMISSION • FLUORESCENCE • FORENSICS • GRATINGS & OEM • OPTICAL SPECTROSCOPY • RAMAN • THIN FILM