



# Analysis of small samples

Various accessories have been developed to accommodate small and odd sample sizes and shapes.

The requirements are : 1) to maintain the vacuum  
2) to assure that the sample at least covers the anode.

For small samples where usually only 1 measurement is possible the HDD mode removes the need for presetting the HV and is an essential tool.



Key points: GD can analyse samples down to 2 mm

**In the USA:**  
Jobin Yvon Inc.  
3880 Park Avenue  
Edison, NJ 08820-3012  
Toll-Free: +1-732-494-8660  
Fax: +1-732-549-5125  
E-mail: info@jobinyvon.com  
www.jobinyvon.com

**In France:**  
Jobin Yvon S.A.S.  
16-18 rue du Canal  
91165 Longjumeau cedex  
Tel: +33 (0) 1 64 54 13 00  
Fax: +33 (0) 1 69 09 07 21  
E-mail: info@jobinyvon.fr  
www.jobinyvon.fr

**In Japan :**  
Horiba Jobin Yvon CO., LTD.  
Higashi-Kanda Daiji Building  
5<sup>th</sup> Floor Higashi-Kanda, Chiyoda-ku, Tokyo  
Tel: +81 (0) 3-3861-8231  
Fax: +81 (0) 3-3861-8253

**Germany:** +49 (0) 89 4623 17-0  
**Italy:** + 39 0 2 57603050  
**UK:** +44 (0) 20 8204 8142

